


⚠️ WARNINGS

- Read, understand and follow all instructions. ⚠️ **Yellow Warning Sign indicates Tensioner is in RELEASE MODE, DO NOT transport occupant when visible.**
- J-Hooks must be attached to a solid wheelchair frame (no spokes or removable components) and Rear J-Hooks at approximate 30-45 degree angle with the floor.
- Shoulder and Lap Belt (we recommend using both) must go across occupant's shoulder and hips, and not be worn twisted or held away from the occupant's body by wheelchair components.



SECURING

(A) REAR

1. Attach shoulder belt pin connector to pin on wall side lap belt, then secure seatbelts onto the storage pin under flip-up seat (A).
2. Backup wheelchair to park zone (B).
3. Pull paddle handle (C) to release rear retractors and activate XPress Timer.
4. Attach both rear J-hooks to rear wheelchair frame (D).
5. Backup wheelchair within park zone (E) and against side bumper or closest proximity.

(B) FRONT

- ⚠️ Make sure Tensioner is in release mode ⚠️ (F).
6. Attach front J-hook to front wall side wheelchair frame (J). For scooters, loop J-hook around front wheel frame, and then insert J-hook (facing out) through O-ring (K).
7. Press button (G), pull down tension handle (H), release button, and crank handle until full tension is reached and clutch slips. Make sure Warning Sign ⚠️ (F) is NOT visible.
8. Lock wheels (L) or power off chair.

(C) SEATBELTS

9. Extend shoulder belt and wall side lap belt across occupant's shoulder and hip (M); then buckle them to the aisle side lap belt (N).
- ⚠️ Buckled seatbelts should be against occupant's hipbone (P).
10. Occupant and wheelchair are secured and ready to go.

RELEASING

(C) SEATBELTS

1. Detach shoulder belt pin connector from pin on lap belt (M).
2. Unbuckle lap belt (N) and place all seatbelts in storage.
- ⚠️ If needed, press seatbelt release button (Q).

(B) FRONT

3. Press button (G), push back tension handle (H) and release button to release front J-hook tension ⚠️ (F).
4. Remove front J-hook from wheelchair frame (J) or scooter (K). Place J-hook in storage ring (R).
5. Unlock wheels (L) or power on chair.

(A) REAR

6. Pull paddle handle (C) to release rear retractors and activate XPress Timer.
7. Move wheelchair forward, remove rear J-hooks from rear wheelchair frame and place J-hooks in storage (D).
8. Occupant and wheelchair are released and ready to go.